


DICOM FILE UPLOAD TO THE BLUEPRINT CLOUD

CONTACT | **SIGN IN** REGISTER

STEP 1: Navigate to shoulderblueprint.com and click on "Sign In" in the upper right corner.





STEP 2: Within the Web Application, navigate to and click on "Dicom Upload" from the left-hand menu.

STEP 3: Select the surgeon you want to receive the files.

A. Once the surgeon is selected, click the green "Confirm" button at the bottom.

Confirm



STEP 4: Drag & Drop on the screen or upload the patient's DICOM file from your computer.

IMPORTANT: Files must be unzipped, extracted or uncompressed to be uploaded. Drag & Drop functionality is **not** compatible while using Internet Explorer.



STEP 7: Once the upload is complete, the files are automatically pre-processed for 3D reconstruction errors and sent to the surgeon's BLUEPRINT account.

Upload Succeeded



Upload a New Case

A. Pending successful 3D reconstruction, the next time the surgeon logs in to their BLUEPRINT software, the case will be "Ready to plan," error-free.

STEP 6: The DICOM files begin to upload into the BLUEPRINT Cloud.

STEP 5: After the DICOM files are selected, click "De-Identify and Upload the Files."

